

Network and Modality Worksheet

This worksheet is used to help you to identify various elements of your network setup.
By identifying the elements of your imaging network you can significantly reduce your setup and implementation time.

Picture Archiving and Capture System Settings:

Server Settings

Host Name	
IP Address	
Application Entity Title This is sometimes referred to as the called application entity title (called AET)	
Primary Modality Work list Port This is sometimes referred as the C-Find Dimse command	
Primary Modality Storage Port This is sometimes referred as the C-Store Dimse command	
Storage Path This is the UNC path for the storage of received dicom objects. Please note that this is a path that your EHR/EMR should be able to access	

Clients and Modalities

Please use this table to identify each modality on your network

Name	IP Address	Client AET	Called AET	MWL Port	Storage Port